AUSTRIA

Bauerfeind Ges.m.b.H.
Hainburger Straße 33
1030 Wien
P +43 (0) 800 4430-130
F +43 (0) 800 4430-131
E info@bauerfeind.at

FRANCE

Bauerfeind France S.A.R.L.
B.P. 59258
95957 Roissy CDG Cedex
P +33 (0) 1 4863-2896
F +33 (0) 1 4863-2963
E info@bauerfeind.fr

SLOVENIA

Bauerfeind d.o.o.
Dolenjska cesta 242 b
1000 Ljubljana
P +386 (0) 1 4272-941
F +386 (0) 1 4272-951
E info@bauerfeind.si

BENELUX

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Polyurethan (PUR),
Polyethylenterephthalat (PET),
Polyester (PES)

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 english

Dear patient,

please read and observe these instructions for use carefully and contact your physician, medical retailer or our Customer Service department at Tel. +49 (0)800 0 01 05 30 (free within

Germany) or at info@bauerfeind.com if you have any questions.

Intended purpose

ErgoPad® weightflex 2 is a medical device for the feet. It is an orthopedic orthosis that is suitable for use in treating and preventing foot complaints. The ErgoPad weightflex 2 encourages the natural gait pattern by providing flexible / central / stable arch support and therefore increases the wearing comfort of your everyday and business shoes.

Indications

The ErgoPad weightflex 2 foot orthoses effectively support treatment in case of foot abnormalities and problems (mild to pronounced arch decreased / spray foot or mild pes cavus / high arched foot).

Risks of using this product

 Caution: before using / applying a foot orthosis, you should consult a physician or specialist*. Only a specialist* may carry out the fitting of the product and provide instruction.

If the available internal volume of the shoe is not sufficient for the foot and the selected foot orthosis, or if the foot orthoses are worn in shoes of an incorrect length and width (or shoes with inadequate lateral guidance in the case of open shoe types), problems relating to poor fit may occur (e.g. a feeling of constriction or pressure, rubbing).

If the feet are positioned on the edges of the foot orthosis, blisters and calluses may form on the soles of the feet. This may also occur if top covers or cushioning material develop holes, creases, or break due to mechanical overloading. You must stop wearing foot orthoses if they are damaged in this way!

Contraindications

People with diabetes and other conditions that affect sensation in the feet and people with rheumatoid foot deformities must seek medical advice before using this product to ensure optimum treatment.

Application instructions

Shoes suitable for wearing with foot orthoses must have sufficient internal volume for both the foot and the selected foot orthosis (adequate shoe length / width and, if necessary, a removable, sufficiently thick cushioned sole, as well as lateral guidance for the foot orthoses in the case of open shoe types). An uncomfortable fit

that constricts the foot must always be avoided. **a**

When selecting the correct size of foot orthosis, the packaging information serves merely as a guide. Just like correctly fitting shoes, foot orthoses must be at least 1.0 cm longer than the foot, when the foot is in a weight-bearing, standing position. Please check the length allowance for both feet, as they may differ in length. In case of doubt, base your size selection on the longer foot. **b**

The foot orthosis must be positioned straight in the shoe and have as much surface contact as possible with the shoe insole. Lateral contact between the foot orthosis and the shoe upper helps reduce the risk of pressure marks or blistering on the skin of the soles of the feet. **c**

Foot orthoses must always be worn as a pair! Any exceptions to this must be stipulated by the treating physician.

Please regularly check your feet for any pressure marks, redness or other changes caused by the foot orthosis. If you have any questions, please contact a specialist or your physician.

Adjustment phase: As with many medical aids, the ErgoPad weightflex 2 may feel somewhat strange to wear at first. We therefore recommend that you start by wearing it for one or two hours per day and slowly increase the wearing period each time until the foot orthosis feels comfortable. Normally, it should not take much longer than a week to adjust to wearing the foot orthosis. If you still find it uncomfortable to wear after this period, please contact a specialist.

Foot orthoses must be checked by a specialist after they have been worn for approximately one year to ensure that they are still functioning properly.

Cleaning instructions

The ErgoPad weightflex 2 can be cleaned by hand using soap and lukewarm water. Then leave the foot orthoses to air dry. Do not use any external heat sources (e.g. heaters, hairdryers or direct sunlight). The foot orthoses cannot be washed in a washing machine.

Mild disinfectants suitable for human skin may be applied sparingly (sprayed on). However, this may lead to discoloration of the top cover materials.

Part of the body this product is used for

The foot orthoses must be worn in suitable footwear and with socks or stockings on your feet.

Notes on reuse

This product is intended to be fitted to one patient only.

Warranty

The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. The product must be cleaned before submitting it in case of warranty claims. Warranty may be limited or excluded if the instructions on how to handle and care for the ErgoPad weightflex 2 have not been observed.

Warranty is excluded if:

- The product was not used according to the indication
 - The instructions given by the specialist were not observed
 - The product was modified arbitrarily
- Discoloration, in particular fading of the top covers of the foot orthoses, may occur during their use. This is dependent on a number of factors, such as moisture levels and sock / stocking materials. This is not covered by warranty.

Disclaimer

Do not diagnose yourself or decide on which medication to take unless you are a specialist*. Before first using our medical device, please seek the advice of a physician or specialist* as this is the only way to assess the effect of our product on your body and determine the potential risks of using this product due to your personal constitution. Follow the advice of this specialist* as well as all instructions in this document or its online version, including any excerpts (as well as texts, images, graphics etc.).

If you still have doubts after your consultation with the specialist*, please contact your physician or distributor, or contact us directly.

Duty to report

Due to regional legal regulations, you are required to immediately report any serious incident involving the use of this medical device to both the manufacturer and the responsible authority. Our contact details can be found on the back of this brochure.

Disposal

Once you have finished using the product, please dispose of it according to local specifications.

Material content

Polyurethane (PUR),
Polyethylene terephthalate (PET),
Polyester (PES)

MD – Medical Device

UDI – DataMatrix Unique Device Identifier UDI

SPECIALISTS*

Modification instructions

Fitting:

For an optimum fit, there should be a length allowance of approximately 1 cm between the foot and the foot orthosis and / or shoe in a standing position. **b**

If the foot orthoses are not the right size for your shoes, you can cut them down or reshape them using thermoforming. Please follow the instructions below.

Cutting:

1. Take the inlays out of your shoes and place them on your ErgoPad weightflex 2.
2. Align the heel sections so that they are positioned precisely on top of one another and draw around the inlay on your ErgoPad weightflex 2.
3. Cut along the outline you have drawn. (Make sure you do not cut the foot orthosis too small. We recommend cutting off less material than you need to at first and trying on the foot orthosis before carrying on cutting).
4. Place the ErgoPad weightflex 2 foot orthoses in your shoes (orthoses for left or right feet are marked as such on the back).

Thermoforming **a**

In exceptional cases and having sought specialist advice, it may be necessary to use thermoforming to adapt the foot orthoses to suit the shape of your feet.

1. After cutting down your ErgoPad weightflex 2, preheat your oven to 90 °C / 200 °F.
2. Put the foot orthoses (with the top covers facing downward on a clean surface) in the oven for approximately 2 minutes.

a Risk of burning! Wear protective gloves when taking the foot orthoses out of the oven.

3. As soon as you have taken the foot orthoses out of the oven, place them immediately in your shoes.
 4. Put on your shoes straight away (make sure you are wearing socks or stocking).
 5. Then stand for 2 minutes with your feet shoulder-width apart and pointing forward.
- Your foot orthoses should now be adapted to the shape of your heel.

Version: 2024-11

*A specialist is any person who is authorized according to the state or institutional regulations as well as other provisions relevant to them for fitting and instruction in the use of orthopedic aids. As the manufacturer, we strongly recommend that this person be trained in orthopedics or possess comparable acquired competencies as well as dispose of orthopedic infrastructure.

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